



DafodilTM

Pericardial Bioprosthesis

Mumbai

Location



Meril is located 100 miles North of Mumbai

Vapi



Meril is 100% owned by Bilakhia Family.

CORE OBJECTIVE

To develop, manufacture & distribute clinically relevant, state-of-the-art & best-in-the-class products to alleviate human suffering & improve quality of life.

Global Med-tech Hub For High Quality, Next Generation Technologies

Robust R&D, Ultra-modern, Large Manufacturing & Commitment to Clinical Research



- Large, hi-tech, mfg. capacity
- Dedicated R&D Complex at Meril Park with more than 150 professionals
- Global clinical program to establish safety and efficacy of products

Wide patient access



Global Market Presence



Creating Long Term Partnership



Strong R&D Focus



Highest Quality Control

Strong Organizational Band Width




Commitment to Training & Education



Large Scale Manufacturing

Presence in 4 Large Business Verticals



Cardiovascular
24 Product Families

- Cardio,
- Peripheral & Neuro vascular
- Structural Heart Disease
- Congestive Heart Failure



Orthopedic
10 Product Families

- Total Knee Replacement
- Total Hip Replacement
- Spine implants
- Trauma implants



Endo- Surgery
30 Product Families

- Sutures
- Staplers
- Adhesives
- Hernia Repair
- Women's health



Diagnostics
45 Product Families

- Analyzers
- Reagents
- Biochemistry
- Hematology
- Immunology
- Rapid testing
- POCT


100+ Technologies



- Largest dedicated training academy for Med-Tech in Asia Pac
- Our Vision: Transform patient care through professional education
- 10,000+ HCPs engaged for various educational programs
- Over 1,000 Interventional Cardiologists from over 20 countries have visited us in last few years
- Dedicated TAVR & Complex PCI (CIT) programs are organized
- Live case transmissions & hands-on simulator models for Ortho, Cardio, Surgery




Over a Decade of Growth and Innovation

Established 2006



- Business in 100+ countries
- 4,000+ employees
- DCGI, CE, USFDA, CFDA, KFDA, ANVISA, INVIMA, Russia MoH, TGA, Canada, PMDA Japan

Vascular



Coronary

- Cobalt Chromium Stent
- PTCA Balloon
- Drug Eluting Stent
- Bioresorbable Scaffold
- Drug Coated Balloon



Peripheral

- Peripheral BMS, DES & BRS
- PTA Balloons
- Nitinol Peripheral BMS & DES
- Drug Coated PTA Balloon

Neurovascular

- Flow modification stents
- Embolisation material
- Micro-catheters
- Accessories

Non-Vascular



Structural Heart

- TAVR
- TMVR
- SVR - Biological / Mechanical
- Valvuloplasty Balloons
- Septal Defect Closure Device

Heart failure

- Intra Aortic Balloon Pump
- Left Ventricular Assist Device


Ablation Therapies

- Renal denervation
- Bronchial thermoplasty

- Product and Features Overview
- Architecture
- Durability
- Implantability
- Hemodynamics
- Dafodil-1 First-In Human trial
- Valve Preparation Procedure
- Sizer set
- Ordering Info
- Technical Specification

Dafodil is a **Bovine Pericardium** tissue based **Tri-leaflet Stented Valve**.

Scientifically modeled and bioengineered design to optimize:



Position-Specific design answers the challenging requirements of the Mitral & Aortic position



Aortic
Bioprosthesis




Mitral
Bioprosthesis

Triple Composite design

Bovine Pericardium Tissue

Polyester knit fabric



Triad Frame Design



Support ring-
Acetal Polymer



Commissures-
PET Film Structure



Wireform-
Elgiloy alloy

Sewing ring

with Commissural markers

(Dafodil Mitral has a Silicone Ring)



Triad Frame Design

Commissure : PET Film structure


Support Ring : Acetal Polymer ring

Wireform : Elgiloy Alloy

- Superior Spring Efficiency
- Fatigue resistant Characteristics
- Corrosion Resistant alloy

Polyester Knit Fabric

- Designed to minimize Para-Valvular leak
- Promotes tissue in-growth





Triple Composite Leaflet Design

Advanced Tissue Thickness Matching Technology

Tissue Deflection Testing to match Size and Shape




Better Hemodynamics and Durability





AntiCa⁺™ is Meril's own Patented Anti-calcification process to mitigate the risk of calcification .

Process




Functions

- Chemically extracts phospholipids from cellular components of pericardial tissue &
- Inactivates the Glutaraldehyde resistant organisms.

Process Outcomes

- Reduces bio-burden
- Renders the tissue resistant to calcification
- Increases Valve's Long term Durability






Aortic Valve

Specially designed Commissure posts ensure :

- Ease of Insertion
- Widen Ostia Clearance
- Compliments more space for better knot tying experience
- Differential Sewing ring with contour design



Mitral Valve

- Sewing ring with additional Silicone ring helps in mitigating Para Valvular leak






Position-Specific Design

Unique Dafodil Mitral Delivery System for Ease of Implant at Mitral position


Surgeons can **Actively Modify** the External Stent Post Diameter (E.S.P.D) to reduce its profile aiding Implantation in the Mitral Annulus.

This is possible by **Locking threads** and **Rotator** unique to mitral valves. When Inserted it creates tension onto the black sutures running over the commissures.

External Stent-Post Diameter



Dafodil Mitral Bioprosthesis Resting Position



Reduced E.S.P.D Profile due to Dafodil Mitral Delivery System



Sewing Ring

- Polyester fabric Sewing ring
 - Soft and Pliable
 - Provides optimal natural fit
- Commissural markers
 - Aid in Proper valve orientation during implantation
 - Helps for Uniform placement of sutures



Sewing Ring - Mitral

- Position-specific design answers the challenging requirements of the Mitral position - additional Silicone ring helps in mitigating Para-Valvular leak




Elgiloy alloy wireform

- Provides Exceptional handling characteristics
- Excellent Radio-Opacity

MRI Compatibility

- MR conditional in a magnetic field of upto 3.0 Tesla




Dafodil Aortic and Mitral Valve Placement



Specifically designed for **Supra-Annular** placement ensuring improved Hemodynamics

Triple-Composite design provides larger **Effective Orifice Area (EOA)**

Leaflet Co-aptation helps minimize Intra-Valvular Regurgitation



Excellent Effective Orifice Area & Leaflet Co-aptation

Purpose

The purpose of Dafodil-1 first-in-human trial was to evaluate clinical safety and performance (including hemodynamic parameters) of the Dafodil pericardial bio-prosthesis in patients who underwent aortic or mitral valve replacement

Method

This prospective, multicenter clinical trial enrolled 60 patients (Aortic: 30 patients; Mitral: 30 patients) from seven sites across India.

- Safety endpoints were early (≤ 30 days) and late (> 30 days) mortality and valve-related morbidity.
- Performance End Points and
- Other End Points

Hiremath et al. *Journal of Cardiothoracic Surgery* (2020) 15:140
<https://doi.org/10.1186/s13019-020-01154-7>

Journal of
Cardiothoracic Surgery

RESEARCH ARTICLE

Open Access

Clinical outcomes and hemodynamic performance of Dafodil[™] aortic and mitral pericardial bioprosthesis: 1-year results from Dafodil-1 first-in-human trial

C. S. Hiremath¹, Anil R. Jain², Anurag Garg³, Nirmal Gupta⁴, Yugal K. Mishra⁵, Zile Singh Meharwal⁶, Nityanand Thakur⁷, Atul A. Maslekar⁸ and Naman Shastri^{9*}



Key safety endpoints

- MACE
- All cause mortality
- Myocardial infarction
- All stroke

Performance endpoints

- NYHA
- Hemodynamic performance
- Device success
- Quality of life
- (SF-12 questionnaire)

Other end points

- Cardiovascular mortality
- Stroke and TIA
- Major and minor bleeding
- Acute kidney injury
- Valve thrombosis
- Structural valve deterioration
- Prosthetic valve endocarditis
- Major paravalvular leak
- Conduction disturbances and arrhythmia
- Non structural valve dysfunction
- Aortic/mitral apparatus damage or dysfunction
- Hemolysis
- Explant
- Trial valve related reoperation



Sr. No.	Principal Investigator	Institute	Enrolled	Aortic	Mitral
1.	Dr. Devi Prasad Shetty	Narayana Institute of Cardiac Sciences, Bangalore	8	6	2
2.	Dr. Nirmal Gupta	Sanjay Gandhi Post Graduate Institute, Lucknow	3	1	2
3.	Dr. Nityanand Thakur	Sassoon Hospital, Pune	1	0	1
4.	Dr. C. S. Hiremath	Sri Sathya Sai Institute of Higher Medical Sciences, Bangalore	22	6	16
5.	Dr. Yugal Mishra*	Fortis Escorts Heart Institute, Delhi	3	2	1
6.	Dr. Anil Jain	SAL Hospital & Medical Institute, Ahmedabad	13	8	5
7.	Dr. Anurag Garg	Dr. D Y Patil Hospital, Pune	10	7	3
Total			60	30	30

*Dr. Yugal Mishra is no more associated with the Fortis Escorts Heart Institute, Delhi. Presently, the study is been conducted at Fortis Escorts Heart Institute, New Delhi with Dr. Z. S. Meherwal as Principal Investigator.

Demographics	AVR	MVR
Age years (mean ± SD)	59.83 ± 8.33	48.57 ± 12.63


Aortic Valve Gender Ratio

Mitral Valve Gender Ratio




Dafodil Aortic and Mitral Pericardial Bioprosthesis
n= 60 (30:30 Aortic:Mitral)

Aortic Valve Size Split



Mitral Valve Size Split




Safety	
Freedom from All-cause Mortality	93.33%
Freedom from Valve related mortality	100%
Freedom from Stroke	100%
Freedom from Valve thrombosis	100%
Freedom from Structural Valve deterioration	100%
Freedom from Permanent Pacemaker Implantation	98.33%
Performance	
Device Success	95%

Serious Adverse Events (SAE) Details:

Sr. No.	Patient ID	Type of SAE
1	Patient 1	Congestive cardiac failure followed by death
2	Patient 2	Death due to cardiopulmonary arrest due to low cardiac output due to biventricular dysfunction
3	Patient 3	Septicemia and multi organ failure leading to death
4	Patient 4	Multi-organ dysfunction syndrome leading to death

All Mortality have been adjudicated as
'Not Related to Study Device' by independent Clinical Event Committee (CEC) of Dafodil-1 FIH Trial

Aortic valve performance (n=30)




Echo analysis performed by Independent Central Echo Laboratory (CBCC Global Research)

~ **72.5%** ↓
Reduction in Mean Pressure gradient*

~ **107%** ↑
Increase in EOA values*

Parameter	Valve Size	Pre-procedure	Discharge	1 Year Data	% Increase at 1 Year
Effective Orifice Area (mm ²)	19 (11)	0.72 ± 0.21	1.59 ± 0.67	1.78 ± 0.28	~147%
	21 (13)	0.68 ± 0.10	1.47 ± 0.55	1.83 ± 0.23	~169%
	23 (3)	1.32 ± NA	1.38 ± 0.11	2.32 ± 0.13	~75.7%
	25 (3)	2.48 ± 1.97	1.60 ± 0.24	1.84 ± 0.11	~-25.8% *
	Overall (30)	0.89 ± 0.70	1.53 ± 0.54	1.85 ± 0.27	~107%

* These patients were a case of dominant Moderate AR with Low Gradients and Higher EOA at Baseline (Pre-Implant)



Echo analysis performed by Independent Central Echo Laboratory (CBCC Global Research)


~ 54.4% ↓
Reduction in Mean Pressure gradient*

~ 91.7% ↑
Increase in EOA values*


*Values compared at Pre-procedure and at 12 months

Parameter	Valve Size	Pre-procedure	Discharge	1 Year Data	% Increase at 1 Year
Effective Orifice Area (mm ²)	25 (5)	0.98 ± 0.13	2.50 ± 0.48	2.37 ± 0.44	~141%
	27 (12)	1.20 ± 0.61	3.10 ± 0.64	2.68 ± 0.55	~123%
	29 (3)	1.93 ± 0.24	3.85 ± 1.38	2.86 ± 0.43	~48%
	31 (10)	1.50 ± 1.13	2.67 ± 0.70	2.78 ± 0.42	~85%
	Overall (30)	1.34 ± 0.77	2.95 ± 0.80	2.67 ± 0.48	~99%

Aortic NYHA class



Mitral NYHA class



NYHA class

■ Class IV

■ Class III

■ Class II

■ Class I

Clinical outcomes of Dafodil-1 study at one-year demonstrate **acceptable preliminary safety and performance** of the Dafodil Pericardial Bioprosthesis implanted at aortic and mitral positions.

Moreover, hemodynamic profile of Dafodil Pericardial Bioprosthesis confirmed favorable performance of the valve at aortic and mitral positions


Objective

To thoroughly rinse the Glutaraldehyde sterilant used as storage solution for the bioprosthesis .


Rinse Procedure

- Place the bioprosthesis in the first bowl of sterile, physiological saline. The saline solution should completely cover the bioprosthesis and holder.
- With the bioprosthesis and holder submerged, slowly agitate (gently swirl) back and forth for a minimum of 1 minute.
- Transfer to 2nd rinsing bowl and gently agitate for a minimum 1 more minute.
- Transfer to 3rd rinsing bowl and gently agitate for a minimum 1 more minute.
- Finally, keep the bioprosthesis and holder in a 4th bowl of saline until implantation to prevent the tissue from drying.

500ml Saline in each Bowl



Aortic sizer set




Mitral sizer set



Available Bioprosthesis Diameter (mm)	Reference Catalogue	
	Aortic Model	Mitral Model
19	DDL19A	N/A
21	DDL21A	N/A
23	DDL23A	DDL23M
25	DDL25A	DDL25M
27	DDL27A	DDL27M
29	N/A	DDL29M
31	N/A	DDL31M


Available Sizer Set	Reference / Catalogue Number
Aortic sizer set	DFA1927
Mitral sizer set	DFM2331

Technical Specification Aortic Valve



Aortic	Size	19	21	23	25	27
A. Tissue Annulus Diameter		20	22	24	26.5	28.5
B. Internal Stent Post Diameter (Base)		18	20	22	24	26
C. External Sewing ring Diameter		27	28	31	33	36

Technical Specification Mitral Valve



Mitral	Size	23	25	27	29	31
A. Stent Diameter (Wire Form)		23	25	27	29	29
B. External Stent Post Diameter (Base)		26.5	28.5	29.5	32.5	32.5
C. External Stent Post Diameter (Tip)		27	29	31	33	33
D. External Sewing Ring Diameter		33	36	38	41	43

Storage Temperature	10° C to 25 ° C (50-77° F)
Storage Solution	Glutaraldehyde
MRI Safety	MR conditional in a magnetic field of upto 3.0 Tesla
Recommendation	Use only Dafodil Aortic/Mitral Sizer Set with Dafodil Aortic/Mitral bioprosthesis

Thank you